



# MF II Fixer

## Premium Universal Fixer

Formulated for consistency and stability with an affordable price

**Versatile Medical Imaging fixer, features high performance with low odor**

### Premium Quality

**FR MF2 Universal Fixer Replenisher** is a two-part high performance fixer concentrate formulated to deliver rapid film drying and compatible with a variety of silver recovery methods. Suitable for RP, extended and manual processing. Complete compatibility with all brands of general diagnostic, laser and mammography imaging films.

#### Technical Information (working Solution)

|                     |                |
|---------------------|----------------|
| pH (at 80°F / 27°C) | 4.20 to 4.30   |
| Specific Gravity    | 1.080 to 1.090 |

### Quality you can count on

Extended chemical stability coupled with consistent batch-to-batch quality, ensures optimal image quality. Low Hypo retention and unexposed silver provides excellent archival characteristics.

### Quick Easy Mixing

Two concentrate parts facilitate fast solution mixing and minimize the chance of mixing errors, either manually or in an automated chemical mixer.

Concentrate bottle templates are available for the popular White Mountain IS199 mixer.



### Low Maintenance

A formulation that resists crystallization results in a clean running solution that minimizes buildup on fixer racks and tanks. The result is reduced maintenance.

### Low Odor

Odor causing acetic acid and sulfur dioxide emissions are significantly reduced, due to an advanced design formulation.

### High Quality Images

**FR MF2 Universal Fixer Replenisher** yields high image quality. It yields a pleasant blue image tone which reduces eye fatigue and improved contrast perception. Fewer artifacts further enhance image quality and reduces costly re-examinations.

### Order in three convenient sizes

|                      |                      |
|----------------------|----------------------|
| To Make 4 Gallons:   | Catalog # 44-602     |
| To Make 10 Gallons:  | Catalog # 44-640     |
| To Make 300 Gallons: | Catalog # 42-798-300 |

### Easy Shipping

Customer shipment transportation surcharges are eliminated due to the non-hazardous formulation designation.

For additional technical information on **FR MF2 Universal Fixer Replenisher** or high quality general medical or mammography developers in the **FR** processing family contact:

**FR Chemical, Inc.**  
300 East Sandford Blvd. 888-634-4646  
Mount Vernon, NY 10550-4626 Fax 914-667-7145