

High Quality Automatic X-Ray Developers

Formulated for optimum image quality, incorporates important health, safety and environmental benefits

Low Odor
Use with all automatic processors
Non-Hazardous & Non-Corrosive

Quick Easy Mixing
Clean working 2-part formulation
Eliminates Hazmat shipping charges



HiC-ES Kodak RP Equivalent Developer & Replenisher

A unique two-part developer concentrate that delivers outstanding diagnostic quality and is ideal for mixed film environments. A multi-purpose performer with all brands of medical imaging films and compatibility with all processing cycles (45, 90 & 180 seconds). Optimum image contrast, coupled with a highly desired blue image tone, results in excellent visualization of subject detail.

Catalogue # 44-505
Catalogue # 44-504 ES
Catalogue # 42-774

Makes 4 Gallons
Makes 10 Gallons
Makes 200 Gallons



ECO Plus Medium Activity Developer & Replenisher

A multi-purpose performer with all brands of medical imaging films, including Agfa, Fuji, Kodak & Konica. The unique two-part developer concentrate is formulated to deliver medium activity while providing consistent diagnostic quality with both general medical and laser imaging films. Complete compatibility with 90 & 180 second processing cycles.

Catalogue # 44-506
Catalogue # 42-776

Makes 10 Gallons
Makes 200 Gallons



MR 2 Mammography Developer & Replenisher

Formulated as a high activity developer, enhanced contrast and optimum speed is achieved with Agfa HDR-C and HDR-C Plus, Fuji AD-M, Kodak MR 2000 and Konica CM-H films.

Ideal for both dedicated and non-dedicated environments. The two-part formulation is non-corrosive and non-hazardous.

Catalogue # 44-508

Makes 10 Gallons

MED-DENT Pre-Mixed Developer & Fixer

Med-Dent solutions have been specifically designed to yield the highest viewbox quality with all brands of Medical and Dental X-Ray films, while providing pre-mixed user convenience. Complete compatibility with all processing cycles and a top performer in self-contained closed loop systems. Acetic acid free developer formula provides reduced odor. A specially formulated fixer formula further reduces chemical odor.

Catalogue # 2500DF	2 each - 1 Gallon Developer bottles 2 each - 1 Gallon Fixer bottles
Catalogue # 2501D	4 each - 1 Gallon Developer bottles
Catalogue # 2502F	4 each - 1 Gallon Fixer bottles



MF II Universal Fixer Replenisher

This premium universal fixer is a two-part high performance fixer concentrate formulated to deliver rapid film drying and minimize emulsion pickoff. It is suitable for RP, extended and manual processing and provides complete compatibility with all brands of general diagnostic, laser and mammography imaging films.

Extended chemical stability coupled with consistent batch to batch quality ensures optimal quality. Low hypo retention provides excellent archival characteristics. Compatible with a wide variety of silver recovery methods.

Catalogue # 44-602	Makes 4 Gallons (2 x 26 oz Bottles / 2 x 6 oz Bottles)
Catalogue # 44-640	Makes 10 Gallons (2 x 1 gal Bottles / 2 x 16 oz Bottles)
Catalogue # 42-798-300	Makes 300 Gallons (1 - 55 gal Drum / 2 x 3.8 gal)



Manual (Hand Tank) Long Life Developer & Replenisher

Ideal for low volume medical and dental imaging. A long tank life is achieved by an exclusive controlled resistance system.

- For use as a developer, one part of concentrate is diluted with four parts of water.
- For use as a replenisher, one part of concentrate is diluted with three parts of water.

The non-corrosive & non-hazardous formulation eliminates transportation surcharges.

Catalogue # 710810	Makes 6 x 1 Gallons (6 x 26 oz Bottles)
Catalogue # 710820	Makes 4 x 5 Gallons (4 x 1 gal Bottles)



Also Available Through FR Chemical

White Mountain Imaging - Film Specific Mammography Developers



T₂ - AG Developer

A high activity, film specific mammography developer, which achieves optimum contrast with Agfa HDR-C and Agfa HDR-C Plus films in both standard and extended processing applications. State of the art design assures the cleanliness and sensitometric stability.

- Also Recommended for Agfa MR-5 II and Agfa HT
- Suitable for standard and extended processing

Catalogue # 4010D-AG

Makes 10 Gallons



T₂ - EX Developer

Specifically designed to provide the highest contrast, without sacrificing detail while maintaining minimum patient exposure with Kodak Min-R 2000 film. Like all T₂ Developers, T₂-EX is non-corrosive and non-hazardous.

- Suitable for standard and extended processing

Catalogue # 4010D-EX

Makes 10 Gallons



T₂ - ADM Developer

Rounds out the family of film specific mammography developers. Formulated to achieve optimum imaging characteristics with Fuji AD-M and Fuji UM-MA mammography film. It is clean working and sensitometrically stable.

- Suitable for standard and extended processing

Catalogue # 4010D-ADM

Makes 10 Gallons

VERSACLEAN Developer Systems Cleaner

A superior cleaner, Versaclean is environmentally safe, easy to use and requires no special procedures. Prior to discharging, check with your local authorities for any additional regulations that may be applicable in your community.

Catalogue # 11-030

Makes up to 12 Gallons

VERSACLEAN Fixer Systems Cleaner

This Non-Chromate concentrate systems cleaner is ready-to-use and environmentally safe.

Catalogue # 11-035

4 x 1 gallon per case
Makes up to 12 gallons

Developer Starter

FR Starter has a neutral pH and may be used with all x-ray developers. For general medical x-ray films use 1.6 oz per gallon. (3.2 oz for Kodak and Agfa Mammography film).

Catalogue # 42-900

6 x 1 quart per case